

## Mobile CR Installation Kit

Retrofits most vehicles



Pictured here in a Chevy HHR



Can be customized to install in the front or back of any vehicle!

### Mobile CR Installation Kit

This unique kit includes a custom designed fiberglass enclosure that creates a desk top surface. Underneath the enclosure is shock absorbing material and tie down straps that protect the CR Reader from moving. The laptop is mounted and locked on the enclosure for added protection. To power the CR a 1500W Pure Sine Wave Inverter is provided along with a battery isolator, circuit breaker and on/off switch.

A detailed wiring diagram and step by step instructions make installing this kit quick and easy!

Also available: **DICOMBridge Software**

**DICOMBridge sends CR Images from Your Vehicle 3-5X's faster!**

#### *Product Features:*

- Increases wireless transmission speed of images
- Routes images to multiple destinations
- Compatible with all CR Systems: Kodak, Fuji, Agfa, etc.
- Creates single or multi-patient CDs & DVDs
- Auto routing capabilities
- Compatible with all major wireless carrier's AirCards